

**SKS**  
Dental Group

**SKS**<sup>®</sup>  
Dental Group



*máxima calidad a precio*  
*en toda nuestra línea de pro*

**SKS**  
Dental Group

**SKS**  
Central  
Spain



# GENERAL CATALOGUE

**SKS**  
Orthodontic

**AS**  
Technology



962 425 366



[www.sksdental.es](http://www.sksdental.es)



[info@sksdental.es](mailto:info@sksdental.es)




# *The perfect ally*

## Quality dental products

**DKS** Dental Group, providing import, export and distribution of dental products. We are a young and professional team.

**DKS** Dental Group exports to Central America, North Africa, East Asia, Northern Europe and our neighbouring country, Portugal.




---

# for the dental professional

## Our Philosophy

Our Philosophy has always been the same: *“be an ally for the dental professional”*. For this reason, we strive every day to become their *“perfect ally”*.





## SKS Orthodontics

---

Premium Kit Boxes.....	4
Metal Brackets .....	6
Aesthetic brackets .....	10
Tubes and Bands .....	13
Archwires .....	18
Elastic bands.....	23
Accessories.....	25
Instruments.....	32
Ligatures and wires.....	38

## As Technology

---

Ultrasonic Bath.....	40
Water distiller.....	41
Gutta Percha Cutter .....	41
Sealing machine.....	42
Contra-angle handpieces and Turbine/Attachment .....	43

## Prime Dental

---

Composites.....	44
Acid etch.....	45
Adhesives / Cements.....	46
Orthodontics.....	47

## Disposable


---

Cotton Rolls / Disposable Swabs.....	48
--------------------------------------	----

## Rike Toys

---

Cuddly toys.....	50
Backpacks .....	54
Toothbrushes .....	57
Pencils .....	58
School pencil cases.....	59
Erasers / Best Patient Certificates.....	60
Stickers / Ink stamps .....	61
Chests and Tooth Holder Pendant.....	62
Tooth-holder Mice / Stress Relief Toys.....	63
Dental Wristbands / Hourglasses.....	64
Slinky toys / Bouncy balls.....	65
Dental Hand Mirrors.....	66



# KITS Premium BOXES

MINT SELF-LIGATING KIT

MINT SELF-LIGATING KIT *Plus*

ICEBERG AESTHETIC SELF-LIGATING KIT

ANDRY WHITE CERAMIC KIT


SAPPHIRE NINA PLUS KIT



## Mint Self-ligating Kit

Ref. ME052MBTKIT  
(MBT 0.022")

Ref. ME052ROTHKIT  
(ROTH 0.022")



## Mint Plus Self-ligating Kit

Ref. ME052MBTKIT-P  
(MBT 0.022")

Ref. ME052ROTHKIT-P  
(ROTH 0.022")





## Each kit includes:

- 10 Cases of brackets
- 1 Xpress Light 3500 mW LED lamp (white or black)
- 1 Straight Pliers TC SKS Instruments
- 1 Distal End Cutter TC SKS Instruments
- 1 pair Weingart Pliers TC SKS Instruments
- 1 Buccal Tube Tweezers SKS Instruments

### Iceberg Aesthetic Self-ligating Kit

Ref. ME052MBT-STKIT  
(MBT 0.022")

Ref. ME052ROTH-STKIT  
(ROTH 0.022")



### Andry White Ceramics Kit

Ref. ME046MBTKIT  
(MBT 0.022")

Ref. ME046KIT  
(ROTH 0.022")

Ref. ME047KIT  
(ROTH 0.018")




### Sapphire Nina Plus Kit

Ref. ME077ZKIT  
(MBT 0.022")

Ref. ME076ZKIT  
(ROTH 0.022")

Ref. ME078ZKIT  
(ROTH 0.018")





## Metal Brackets



### Celu Series Mini Bracket

- Low-friction mini metal bracket.
- Manufactured with the latest MIM (metal injection moulding) technology.
- Ultra-fine design with a low profile, reducing bracket thickness and enhancing patient comfort.
- Faster tooth sliding, with less friction (shorter treatments).
- Vertical marking for easier positioning.
- Integrated MESH-calibre base with surface sandblasting treatment for excellent bonding.
- Rounded edges.
- Hooks on 3, 4 and 5.

#### Kit containing 20 pcs

Prescription	Ref.
.018 ROTH	ME043
.022 ROTH	ME042
.022 MBT	ME042MBT

ROTH prescription										
	Upper/Right					Upper/Left				
Size	3.0 mm	3.0 mm	3.4 mm	3.0 mm	3.7 mm	3.7 mm	3.0 mm	3.4 mm	3.0 mm	3.0 mm
Torque	-7	-7	0	8	12	12	8	0	-7	-7
Ang	0	0	11	9	5	5	9	11	0	0
in/out	0.6 mm	0.6 mm	0.6 mm	1.0 mm	1.0 mm	1.0 mm	1.0 mm	0.6 mm	0.6 mm	0.6 mm
	5	4	3	2	1	1	2	3	4	5
	●	●	●	●	●	●	●	●	●	●
	●	●	●	○	○	○	○	●	●	●
	5	4	3	2	1	1	2	3	4	5
Size	3.2 mm	3.2 mm	3.4 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.4 mm	3.2 mm	3.2 mm
Torque	-22	-17	-11	0	0	0	0	-11	-17	-22
Ang	0	0	5	0	0	0	0	5	0	0
in/out	0.8 mm	0.8 mm	0.8 mm	1.09 mm	1.09 mm	1.09 mm	1.09 mm	0.8 mm	0.8 mm	0.8 mm
	Lower/Right					Lower/Left				

MBT prescription										
	Upper/Right					Upper/Left				
Size	3.0 mm	3.0 mm	3.4 mm	3.0 mm	3.7 mm	3.7 mm	3.0 mm	3.4 mm	3.0 mm	3.0 mm
Torque	-7	-7	0	10	17	17	10	0	-7	-7
Ang	0	0	8	8	4	4	8	8	0	0
in/out	0.6 mm	0.6 mm	0.6 mm	1.0 mm	0.6 mm	0.6 mm	1.0 mm	0.6 mm	0.6 mm	0.6 mm
	5	4	3	2	1	1	2	3	4	5
	●	●	●	●	●	●	●	●	●	●
	●	●	●	○	○	○	○	●	●	●
	5	4	3	2	1	1	2	3	4	5
Size	3.2 mm	3.2 mm	3.4 mm	2.4 mm	2.4 mm	2.4 mm	2.4 mm	3.4 mm	3.2 mm	3.2 mm
Torque	-17	-12	-6	-6	-6	-6	-6	-6	-12	-17
Ang	0	0	3	0	0	0	0	3	0	0
in/out	0.5 mm	0.5 mm	0.6 mm	1.0 mm	1.0 mm	1.0 mm	1.0 mm	0.6 mm	0.6 mm	0.6 mm
	Lower/Right					Lower/Left				

#### Replacements 5 pcs





## Mint Self-Ligating Bracket



- Low-friction self-ligating bracket.
- Adopts “passive slide” to provide a constant force.
- Faster tooth sliding, with less friction (shorter treatments).
- Removable hook.
- Soft rounded edges.
- Sliding clip.
- Laser- printed tooth numbering.
- Passive system.

### Set of 28 pcs Tubes included


Prescription	Ref.
.022 ROTH	ME052ROTH
.022 MBT	ME052MBT

ROTH 022 prescription			
	TOOTH	TORQUE	ANG
MAXILLARY	Central	+12°	+5°
	Lateral	+8°	+9°
	Cuspid with Hook	-2°	+13°
	1 <sup>st</sup> Bicuspid with Hook	-7°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-7°	0°
MANDIBULAR	Anterior	-1°	0°
	Cuspid with Hook	-11°	+7°
	1 <sup>st</sup> Bicuspid with Hook	-17°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-22°	0°

### Replacements 1 pc



MBT 022 prescription			
	TOOTH	TORQUE	ANG
MAXILLARY	Central	+17°	+4°
	Lateral	+10°	+8°
	Cuspid with Hook	-7°	+8°
	1 <sup>st</sup> Bicuspid with Hook	-7°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-7°	0°
MANDIBULAR	Anterior	-6°	0°
	Cuspid with Hook	-6°	+3°
	1 <sup>st</sup> Bicuspid with Hook	-12°	+2°
	2 <sup>nd</sup> Bicuspid with Hook	-17°	+2°



## Metal Brackets



### Mint Self-Ligating Bracket Plus

Low-friction self-ligating bracket. Adopts “passive slide” to provide a constant force.

Faster tooth sliding, with less friction (shorter treatments).

Ultra-fine design with a low profile, reducing bracket thickness and enhancing patient comfort.

- Soft rounded edges.
- Sliding clip.
- Passive system.
- Laser-printed tooth numbering.
- Ultra-fine low-profile design.

#### Set of 28 pcs Tubes included

Prescription	Ref.
.022 ROTH	ME052ROTH-P
.022 MBT	ME052MBT-P

ROTH 022 prescription			
	TOOTH	TORQUE	ANG
MAXILLARY	Central	+12°	+5°
	Lateral	+8°	+9°
	Cuspid with Hook	-2°	+13°
	1 <sup>st</sup> Bicuspid with Hook	-7°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-7°	0°
MANDIBULAR	Anterior	-1°	0°
	Cuspid with Hook	-11°	+7°
	1 <sup>st</sup> Bicuspid with Hook	-17°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-22°	0°

#### Replacements 1 pc



MBT 022 prescription			
	TOOTH	TORQUE	ANG
MAXILLARY	Central	+17°	+4°
	Lateral	+10°	+8°
	Cuspid with Hook	-7°	+8°
	1 <sup>st</sup> Bicuspid with Hook	-7°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-7°	0°
MANDIBULAR	Anterior	-6°	0°
	Cuspid with Hook	-6°	+3°
	1 <sup>st</sup> Bicuspid with Hook	-12°	+2°
	2 <sup>nd</sup> Bicuspid with Hook	-17°	+2°



## Mint Self-Ligating Bracket Active

Interactive self-ligating brackets provide very low friction in the initial stages of treatment and once the archwires have been formed, the interactive clip presses the archwire right into the bracket slot, thus achieving correct dental positioning.

- Manufactured with MIM technology.
- 17-4 PH Stainless steel metal bracket.
- MIM technology, two-part construction, with laser soldering, without flux.
- Extra-slim design, with rounded edges for maximum patient comfort.
- Active locking clip that can be opened and closed with just a probe.
- 80 calibre MESH with laser numbering.

### Set of 28 pcs Tubes included


Prescription	Ref.
.022 ROTH	ME052ROTH-ACT
.022 MBT	ME052MBT-ACT

ROTH 022 prescription			
	TOOTH	TORQUE	ANG
MAXILLARY	Central	+12°	+5°
	Lateral	+8°	+9°
	Cuspid with Hook	-2°	+13°
	1 <sup>st</sup> Bicuspid with Hook	-7°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-7°	0°
MANDIBULAR	Anterior	-1°	0°
	Cuspid with Hook	-11°	+7°
	1 <sup>st</sup> Bicuspid with Hook	-17°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-22°	0°

### Replacements 1 pc



MBT 022 prescription			
	TOOTH	TORQUE	ANG
MAXILLARY	Central	+17°	+4°
	Lateral	+10°	+8°
	Cuspid with Hook	-7°	+8°
	1 <sup>st</sup> Bicuspid with Hook	-7°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-7°	0°
MANDIBULAR	Anterior	-6°	0°
	Cuspid with Hook	-6°	+3°
	1 <sup>st</sup> Bicuspid with Hook	-12°	+2°
	2 <sup>nd</sup> Bicuspid with Hook	-17°	+2°



## Aesthetic brackets



## Andry White Ceramic Bracket Plus

Aesthetic ceramic bracket, made with polycrystalline alumina powders, for a more translucent appearance. CIM technology, manufactured in one piece.

Form rounded polished edges with a diamond disc. Smooth surface, rounded edges, for a more comfortable support for the patient.

- The slot grooves are bevelled for guaranteed precision even in the deepest part of the groove.
- Four horizontal bar base for high tooth retention.
- Hooks on 3, 4 and 5.

### Kit containing 20 pcs

Prescription	Ref.
.022 ROTH	ME046
.018 ROTH	ME047
.022 MBT	ME046MBT

ROTH prescription				
	TOOTH		TORQUE	ANG
MAXILLARY	1   1	Anterior	+12°	+5°
	2   2	2 <sup>nd</sup> Incisor	+8°	+9°
	3   3	Cuspid with Hook	-2°	+10°
	54   45	Bicuspid with Hook	-7°	0°
	54   45	Bicuspid	-7°	0°
MANDIBULAR	21   12	Lower Anterior	-1°	0°
	3   3	Cuspid with Hook	-11°	+7°
	4   4	1 <sup>st</sup> Bicuspid with Hook	-17°	0°
	4   4	1 <sup>st</sup> Bicuspid	-17°	0°
	5   5	2 <sup>nd</sup> Bicuspid with Hook	-22°	0°
	5   5	2 <sup>nd</sup> Bicuspid	-22°	0°

### Replacements 5 pcs



MBT prescription				
	TOOTH		TORQUE	ANG
MAXILLARY	1   1	Anterior	+17°	+4°
	2   2	2 <sup>nd</sup> Incisor	+10°	+8°
	3   3	Cuspid with Hook	-7°	+8°
	54   45	Bicuspid with Hook	-7°	0°
	54   45	Bicuspid	-7°	0°
MANDIBULAR	21   12	Lower Anterior	-6°	0°
	3   3	Cuspid with Hook	-6°	+3°
	4   4	1 <sup>st</sup> Bicuspid with Hook	-12°	+2°
	4   4	1 <sup>st</sup> Bicuspid	-12°	+2°
	5   5	2 <sup>nd</sup> Bicuspid with Hook	-17°	+2°
	5   5	2 <sup>nd</sup> Bicuspid	-17°	+2°

## Nina Sapphire Bracket Plus



Monocrystalline sapphire bracket, with 97% transparency, mimicking the colour of the teeth to achieve the aesthetic desired by the patient.

- CIM technology, manufactured as one piece.
- Contoured design, crystal- polished surface, low friction, achieving optimum sliding and improved performance.
- Robustness and rigidity for precision placement for tooth torque and rotation.
- With indicator (coloured plastic cross) for perfect positioning of brackets.
- Hooks on 3, 4 and 5.

### Kit containing 20 pcs


Prescription	Ref.
.018 ROTH	ME077Z
.022 ROTH	ME076Z
.022 MBT	ME078Z

ROTH prescription				
TOOTH			TORQUE	ANG
MAXILLARY	1   1	Anterior	+12°	+5°
	2   2	2 <sup>nd</sup> Incisor	+8°	+9°
	3   3	Cuspid with Hook	-2°	+10°
	54   45	Bicuspid with Hook	-7°	0°
	54   45	Bicuspid	-7°	0°
MANDIBULAR	21   12	Lower Anterior	-6°	0°
	3   3	Cuspid with Hook	-6°	+3°
	4   4	1 <sup>st</sup> Bicuspid with Hook	-12°	+2°
	4   4	1 <sup>st</sup> Bicuspid	-12°	+2°
	5   5	2 <sup>nd</sup> Bicuspid with Hook	-17°	+2°
	5   5	2 <sup>nd</sup> Bicuspid	-17°	+2°

### Replacements 5 pcs



MBT prescription				
TOOTH			TORQUE	ANG
MAXILLARY	1   1	Anterior	+17°	+4°
	2   2	2 <sup>nd</sup> Incisor	+10°	+8°
	3   3	Cuspid with Hook	-7°	+8°
	54   45	Bicuspid with Hook	-7°	0°
	54   45	Bicuspid	-7°	0°
MANDIBULAR	21   12	Lower Anterior	-6°	0°
	3   3	Cuspid with Hook	-6°	+3°
	4   4	1 <sup>st</sup> Bicuspid with Hook	-12°	+2°
	4   4	1 <sup>st</sup> Bicuspid	-12°	+2°
	5   5	2 <sup>nd</sup> Bicuspid with Hook	-17°	+2°
	5   5	2 <sup>nd</sup> Bicuspid	-17°	+2°



## Aesthetic brackets



## Passive Self-Ligating Aesthetic Bracket **ICEBERG**

Highly transparent bioceramic bracket without discolouration. Ceramic with highest purity aluminium oxide (Al<sub>2</sub>O<sub>3</sub>), with extremely high transparency. The high-purity aluminium effectively reduces discolouration.

- CIM (ceramic injection moulding) technology.
- Low profile, rounded edges, double-polished surface.
- State-of-the-art passive self-ligating technology achieves light force movements that reduce chair time efficiently.
- Central hook in upper canines and premolars.
- Metal bracket in 34-35-44-45, with distal hook.

### BRACKET SLOT MEASUREMENTS

Slot	Width	Depth
0.022"	0.56 mm	0.76mm

### Kit containing 20 pcs

Prescription	Ref.
.022 ROTH	ME052ROTH-ST
.022 MBT	ME052MBT-ST


ROTH 022 prescription			
	TOOTH	TORQUE	ANG
MAXILLARY	Central	+12°	+5°
	Lateral	+8°	+9°
	Cuspid with Hook	-2°	+11°
	1 <sup>st</sup> Bicuspid with Hook	-7°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-7°	0°
MANDIBULAR	Anterior	0°	0°
	Cuspid with Hook	-11°	+7°
	1 <sup>st</sup> Bicuspid with Hook	-17°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-22°	0°

### Replacements 1 pc



MBT 022 prescription			
	TOOTH	TORQUE	ANG
MAXILLARY	Central	+17°	+4°
	Lateral	+10°	+8°
	Cuspid with Hook	-7°	+8°
	1 <sup>st</sup> Bicuspid with Hook	-7°	0°
	2 <sup>nd</sup> Bicuspid with Hook	-7°	0°
MANDIBULAR	Anterior	-6°	0°
	Cuspid with Hook	-6°	+3°
	1 <sup>st</sup> Bicuspid with Hook	-12°	+2°
	2 <sup>nd</sup> Bicuspid with Hook	-17°	+2°





## Non-Convertible Buccal Tube

Easy to place tube with direct cementing. Manufactured with MIM technology for increased clinical efficiency without reducing patient comfort.

- Rounded surfaces.
- With colour marking for identification.
- Reduced occlusal interface and less gum irritation.
- With marking for better orientation while placing, as well as a broad opening.
- Manufactured with 80-micron mesh.



**\*1 Bag of 10.**  
(same quadrant)



**\*2 Box of 40.**  
(10 tubes for each quadrant)

ROTH 0.18 (10 pcs.) *1				
Non-Conv. Ref.	Tooth	Torque	Offset	
ME142-10SD	16	-10	8°	1st Molar
ME142-10SI	26	-10	8°	
ME142-10II	36	-25	5°	
ME142-10ID	46	-25	5°	
ME136-10SD	17	-10	8°	2nd Molar
ME136-10SI	27	-10	8°	
ME136-10II	37	-25	5°	
ME136-10ID	47	-25	5°	


ROTH 0.18 (40 pcs.) *2				
Non-Conv. Ref.	Tooth	Torque	Offset	
ME142-40	6			
ME136-40	7			

ROTH 0.22 (10 pcs.) *1				
Non-Conv. Ref.	Tooth	Torque	Offset	
ME141-10SD	16	-10	8°	1st Molar
ME141-10SI	26	-10	8°	
ME141-10II	36	-25	5°	
ME141-10ID	46	-25	5°	
ME135-10SD	17	-10	8°	2nd Molar
ME135-10SI	27	-10	8°	
ME135-10II	37	-25	5°	
ME135-10ID	47	-25	5°	

ROTH 0.22 (40 pcs.) *2				
Non-Conv. Ref.	Tooth	Torque	Offset	
ME141-40	6			
ME135-40	7			

MBT 0.22 (10 pcs.) *1				
Non-Conv. Ref.	Tooth	Torque	Offset	
ME141-10MBTSD	16	-14	10°	1st Molar
ME141-10MBTSI	26	-14	10°	
ME141-10MBTII	36	-20	0°	
ME141-ID 10MBTID	46	-20	0°	
ME135-10MBTSD	17	-14	10°	2nd Molar
ME135-10MBTSI	27	-14	10°	
ME135-10MBTII	37	-20	0°	
ME135-ID 10MBTID	47	-20	0°	

MBT 0.22 (40 pcs.) *2				
Non-Conv. Ref.	Tooth	Torque	Offset	
ME141-40MBT	6			
ME135-40MBT	7			



## Cementing Tubes



Double convertible tubes available **Request information**

### Convertible Buccal Tube



The tube has been carefully designed by studying the morphology of the vestibular side of the patient in depth, resulting in a small tube that avoids discomfort, occlusal interference and decementation, making it highly convenient for patients.

- Easy to place with direct cementing.
- Colour marking for identification.
- Rounded surfaces.
- Area to accommodate instruments.
- Mark for improved orientation.
- 80-micron mesh.
- Broad opening.

**\*1 Bag of 10.**  
(same quadrant)

ROTH 0.18 (10 pcs.) *1				
Conv. Ref.	Tooth	Torque	Offset	
ME142-C10SD	16	-10	8°	
ME142-C10SI	26	-10	8°	
ME142-C10II	36	-25	5°	
ME142-C10ID	46	-25	5°	

1st  
Molar

ROTH 0.22 (10 pcs.) *1				
Conv. Ref.	Tooth	Torque	Offset	
ME141-C10SD	16	-10	8°	
ME141-C10SI	26	-10	8°	
ME141-C10II	36	-25	5°	
ME141-C10ID	46	-25	5°	

1st  
Molar

MBT 0.22 (10 pcs.) *1				
Conv. Ref.	Tooth	Torque	Offset	
ME141-C10MBTSD	16	-14	10°	
ME141-C10MBTSI	26	-14	10°	
ME141-C10MBTII	36	-20	0°	
ME141-C10MBTID	46	-20	0°	


1st  
Molar

### Band sorter box



- Stackable modular organiser.
- Accurately locate the parts you need at all times
- Made of durable plastic material
- Includes four trays

**Pack of 1**  
Ref. M1-BOX



## Bands

### Anatomical Band with Double Tube and Lower Wing Nut




- All parts with laser pre-welding.
- Optimum anatomical form
- The perfect balance between resilience for adaptation and ductility
- Roughened interior.
- Permanent laser marking.

**Bag of 2**

Band No.	Orthodontics Tooth No.	
	Ref. No. 36	Ref. No. 46
29	M1-29-TD-LR-CP	M1-29-TD-LL-CP
29+	M1-29+-TD-LR-CP	M1-29+-TD-LL-CP
30	M1-30-TD-LR-CP	M1-30-TD-LL-CP
30+	M1-30+-TD-LR-CP	M1-30+-TD-LL-CP
31	M1-31-TD-LR-CP	M1-31-TD-LL-CP
31+	M1-31+-TD-LR-CP	M1-31+-TD-LL-CP
32	M1-32-TD-LR-CP	M1-32-TD-LL-CP
32+	M1-32+-TD-LR-CP	M1-32+-TD-LL-CP
33	M1-33-TD-LR-CP	M1-33-TD-LL-CP
33+	M1-33+-TD-LR-CP	M1-33+-TD-LL-CP
34	M1-34-TD-LR-CP	M1-34-TD-LL-CP
34+	M1-34+-TD-LR-CP	M1-34+-TD-LL-CP
35	M1-35-TD-LR-CP	M1-35-TD-LL-CP
35+	M1-35+-TD-LR-CP	M1-35+-TD-LL-CP
36	M1-36-TD-LR-CP	M1-36-TD-LL-CP
36+	M1-36+-TD-LR-CP	M1-36+-TD-LL-CP

Band No.	Orthodontics Tooth No.	
	Ref. No. 36	Ref. No. 46
37	M1-37-TD-LR-CP	M1-37-TD-LL-CP
37+	M1-37+-TD-LR-CP	M1-37+-TD-LL-CP
38	M1-38-TD-LR-CP	M1-38-TD-LL-CP
38+	M1-38+-TD-LR-CP	M1-38+-TD-LL-CP
39	M1-39-TD-LR-CP	M1-39-TD-LL-CP
39+	M1-39+-TD-LR-CP	M1-39+-TD-LL-CP
40	M1-40-TD-LR-CP	M1-40-TD-LL-CP
40+	M1-40+-TD-LR-CP	M1-40+-TD-LL-CP
41	M1-41-TD-LR-CP	M1-41-TD-LL-CP
41+	M1-41+-TD-LR-CP	M1-41+-TD-LL-CP
42	M1-42-TD-LR-CP	M1-42-TD-LL-CP
42+	M1-42+-TD-LR-CP	M1-42+-TD-LL-CP
43	M1-43-TD-LR-CP	M1-43-TD-LL-CP
43+	M1-43+-TD-LR-CP	M1-43+-TD-LL-CP
44	M1-44-TD-LR-CP	M1-44-TD-LL-CP
44+	M1-44+-TD-LR-CP	M1-44+-TD-LL-CP





## Anatomical Band with Triple Tube and Upper Wing



- All parts with laser pre-welding.
- Optimum anatomical form
- The perfect balance between resilience for adaptation and ductility
- Roughened interior.
- Permanent laser marking.

**Bag of 2**

Band No.	Orthodontics Tooth No.	
	Ref. No. 16	Ref. No. 26
29	M1-29-TT-UR-CC	M1-29-TT-UL-CC
29+	M1-29+-TT-UR-CC	M1-29+-TT-UL-CC
30	M1-30-TT-UR-CC	M1-30-TT-UL-CC
30+	M1-30+-TT-UR-CC	M1-30+-TT-UL-CC
31	M1-31-TT-UR-CC	M1-31-TT-UL-CC
31+	M1-31+-TT-UR-CC	M1-31+-TT-UL-CC
32	M1-32-TT-UR-CC	M1-32-TT-UL-CC
32+	M1-32+-TT-UR-CC	M1-32+-TT-UL-CC
33	M1-33-TT-UR-CC	M1-33-TT-UL-CC
33+	M1-33+-TT-UR-CC	M1-33+-TT-UL-CC
34	M1-34-TT-UR-CC	M1-34-TT-UL-CC
34+	M1-34+-TT-UR-CC	M1-34+-TT-UL-CC
35	M1-35-TT-UR-CC	M1-35-TT-UL-CC
35+	M1-35+-TT-UR-CC	M1-35+-TT-UL-CC
36	M1-36-TT-UR-CC	M1-36-TT-UL-CC
36+	M1-36+-TT-UR-CC	M1-36+-TT-UL-CC

Band No.	Orthodontics Tooth No.	
	Ref. No. 16	Ref. No. 26
37	M1-37-TT-UR-CC	M1-37-TT-UL-CC
37+	M1-37+-TT-UR-CC	M1-37+-TT-UL-CC
38	M1-38-TT-UR-CC	M1-38-TT-UL-CC
38+	M1-38+-TT-UR-CC	M1-38+-TT-UL-CC
39	M1-39-TT-UR-CC	M1-39-TT-UL-CC
39+	M1-39+-TT-UR-CC	M1-39+-TT-UL-CC
40	M1-40-TT-UR-CC	M1-40-TT-UL-CC
40+	M1-40+-TT-UR-CC	M1-40+-TT-UL-CC
41	M1-41-TT-UR-CC	M1-41-TT-UL-CC
41+	M1-41+-TT-UR-CC	M1-41+-TT-UL-CC
42	M1-42-TT-UR-CC	M1-42-TT-UL-CC
42+	M1-42+-TT-UR-CC	M1-42+-TT-UL-CC
43	M1-43-TT-UR-CC	M1-43-TT-UL-CC
43+	M1-43+-TT-UR-CC	M1-43+-TT-UL-CC
44	M1-44-TT-UR-CC	M1-44-TT-UL-CC
44+	M1-44+-TT-UR-CC	M1-44+-TT-UL-CC

## Anatomical Band with Triple Tube and Upper Wing




- All parts with laser pre-welding.
- Optimum anatomical form
- The perfect balance between resilience for adaptation and ductility
- Roughened interior.
- Permanent laser marking.

**Bag of 2**

Band No.	Orthodontics Tooth No.	
	Ref. No. 16	Ref. No. 26
29	M1-29-TT-UR-CP	M1-29-TT-UL-CP
29+	M1-29+-TT-UR-CP	M1-29+-TT-UL-CP
30	M1-30-TT-UR-CP	M1-30-TT-UL-CP
30+	M1-30+-TT-UR-CP	M1-30+-TT-UL-CP
31	M1-31-TT-UR-CP	M1-31-TT-UL-CP
31+	M1-31+-TT-UR-CP	M1-31+-TT-UL-CP
32	M1-32-TT-UR-CP	M1-32-TT-UL-CP
32+	M1-32+-TT-UR-CP	M1-32+-TT-UL-CP
33	M1-33-TT-UR-CP	M1-33-TT-UL-CP
33+	M1-33+-TT-UR-CP	M1-33+-TT-UL-CP
34	M1-34-TT-UR-CP	M1-34-TT-UL-CP
34+	M1-34+-TT-UR-CP	M1-34+-TT-UL-CP
35	M1-35-TT-UR-CP	M1-35-TT-UL-CP
35+	M1-35+-TT-UR-CP	M1-35+-TT-UL-CP
36	M1-36-TT-UR-CP	M1-36-TT-UL-CP
36+	M1-36+-TT-UR-CP	M1-36+-TT-UL-CP


Band No.	Orthodontics Tooth No.	
	Ref. No. 16	Ref. No. 26
37	M1-37-TT-UR-CP	M1-37-TT-UL-CP
37+	M1-37+-TT-UR-CP	M1-37+-TT-UL-CP
38	M1-38-TT-UR-CP	M1-38-TT-UL-CP
38+	M1-38+-TT-UR-CP	M1-38+-TT-UL-CP
39	M1-39-TT-UR-CP	M1-39-TT-UL-CP
39+	M1-39+-TT-UR-CP	M1-39+-TT-UL-CP
40	M1-40-TT-UR-CP	M1-40-TT-UL-CP
40+	M1-40+-TT-UR-CP	M1-40+-TT-UL-CP
41	M1-41-TT-UR-CP	M1-41-TT-UL-CP
41+	M1-41+-TT-UR-CP	M1-41+-TT-UL-CP
42	M1-42-TT-UR-CP	M1-42-TT-UL-CP
42+	M1-42+-TT-UR-CP	M1-42+-TT-UL-CP
43	M1-43-TT-UR-CP	M1-43-TT-UL-CP
43+	M1-43+-TT-UR-CP	M1-43+-TT-UL-CP
44	M1-44-TT-UR-CP	M1-44-TT-UL-CP
44+	M1-44+-TT-UR-CP	M1-44+-TT-UL-CP





## Preformed Archwires

### Niti Archwire Super Elastic Ovoid




Nickel Titanium (NiTi) Wire. Verified arch shape with the strictest quality control.

- Ideal robustness and elasticity.
- Dual mechanical/electrochemical surface polishing to ensure maximum sliding.
- Available in all sizes for both upper and lower arches.
- Ovoid shape

**Pack: 10 pcs**

Ovoid		
Round	Ref. Upper	Ref. Lower
.012"	SE067U	SE067L
.014"	SE068U	SE068L
.016"	SE069U	SE069L
.018"	SE070U	SE070L
.020"	SE071U	SE071L

Ovoid		
Rectangular	Ref. Upper	Ref. Lower
.016"X .016"	SE072U	SE072L
.016"X .022"	SE073U	SE073L
.017"X .022"	SE074U	SE074L
.017"X .025"	SE076U	SE076L
.018"X .022"	SE075U	SE075L
.018"X .025"	SE077U	SE077L
.019"X .025"	SE078U	SE078L
.021"X .025"	SE079U	SE079L



### Niti Archwire Super Elastic Natural

Nickel Titanium (NiTi) Wire. Verified arch shape with the strictest quality control.

- Ideal robustness and elasticity.
- Dual mechanical/electrochemical surface polishing to ensure maximum sliding.
- Available in all sizes for both upper and lower arches.
- Natural Shape

**Pack: 10 pcs**

Natural		
Round	Ref. Upper	Ref. Lower
.012"	SQ150U	SQ150L
.014"	SQ151U	SQ151L
.016"	SQ152U	SQ152L
.018"	SQ153U	SQ153L
.020"	SQ154U	SQ154L

Natural		
Rectangular	Ref. Upper	Ref. Lower
.014"X .025"	SQ163U	SQ163L
.016"X .016"	SQ155U	SQ155L
.016"X .022"	SQ156U	SQ156L
.017"X .022"	SQ157U	SQ157L
.017"X .025"	SQ158U	SQ158L
.018"X .022"	SQ159U	SQ159L
.018"X .025"	SQ160U	SQ160L
.019"X .025"	SQ161U	SQ161L
.021"X .025"	SQ162U	SQ162L



## Niti Archwire Thermal Activate



Dual mechanical and electrochemical surface polishing of archwire to ensure maximum sliding.

Has the properties and performance of the Niti super elastic archwire.

Provides greater elasticity with less stiffness.

Material fully malleable at temperatures below 5°C.

Its thermal properties achieve optimum activation at 37°C. This unique characteristic allows easy placement of the archwire in any type of malocclusion.

**Pack: 5 pcs**

Round	Ovoid	
	Ref. Upper	Ref. Lower
.012"	TE054U	TE054L
.014"	TE055U	TE055L
.016"	TE056U	TE056L
.018"	TE057U	TE057L
.020"	TE058U	TE058L

Rectangular	Ovoid	
	Ref. Upper	Ref. Lower
.016"X .016"	TE059U	TE059L
.016"X .022"	TE060U	TE060L
.017"X .022"	TE061U	TE061L
.017"X .025"	TE063U	TE063L
.018"X .022"	TE062U	TE062L
.018"X .025"	TE064U	TE064L
.019"X .025"	TE065U	TE065L
.021"X .025"	TE066U	TE066L

## Niti Archwire Super Elastic Aesthetic Ovoid




• Properties and performance of the highest quality Nickel Titanium Super Elastic.

• Dual mechanical and electrochemical polishing of archwire to ensure maximum sliding.

**Pack: 5 pcs**


Round	Ovoid	
	Ref. Upper	Ref. Lower
.012"	ST012U	ST012L
.014"	ST013U	ST013L
.016"	ST014U	ST014L
.018"	ST015U	ST015L
.020"	ST016U	ST016L

Rectangular	Ovoid	
	Ref. Upper	Ref. Lower
.016"X .016"	ST017U	ST017L
.016"X .022"	ST018U	ST018L
.017"X .022"	ST019U	ST019L
.017"X .025"	ST021U	ST021L
.018"X .022"	ST020U	ST020L
.018"X .025"	ST022U	ST022L
.019"X .025"	ST023U	ST023L



## Preformed Archwires

### Natural Stainless Steel Archwire




- Heat-treated steel archwire.
- Robustness and elasticity.
- Dual mechanical/electrochemical polishing.
- Available in all sizes.

**Pack: 10 pcs**

Natural		
Round	Ref. Upper	Ref. Lower
.012"	ACN034U	ACN034L
.014"	ACN035U	ACN035L
.016"	ACN036U	ACN036L
.018"	ACN037U	ACN037L
.020"	ACN038U	ACN038L

Natural		
Rectangular	Ref. Upper	Ref. Lower
.016" X .016"	ACN039U	ACN039L
.016" X .022"	ACN040U	ACN040L
.017" X .022"	ACN045U	ACN045L
.017" X .025"	ACN041U	ACN041L
.018" X .022"	ACN046U	ACN046L
.018" X .025"	ACN042U	ACN042L
.019" X .025"	ACN043U	ACN043L
.021" X .025"	ACN044U	ACN044L

### Ovoid Stainless Steel Archwire



- Heat-treated steel archwire.
- Robustness and elasticity.

**Pack: 10 pcs**

Ovoid		
Round	Ref. Upper	Ref. Lower
.012"	AC034U	AC034L
.014"	AC035U	AC035L
.016"	AC036U	AC036L
.018"	AC037U	AC037L
.020"	AC038U	AC038L

Ovoid		
Rectangular	Ref. Upper	Ref. Lower
.016" X .016"	AC039U	AC039L
.016" X .022"	AC040U	AC040L
.017" X .022"	AC045U	AC045L
.017" X .025"	AC041U	AC041L
.018" X .022"	AC046U	AC046L
.018" X .025"	AC042U	AC042L
.019" X .025"	AC043U	AC043L
.021" X .025"	AC045U	AC045L

## Reverse Curve Niti Archwire



- High-quality Nickel-Titanium alloy.
- Progressive increase of residual forces, full recovery.
- Pre-compensated Anti-Spee shape.
- Easy bending with lower load.

**Pack: 2 pcs**

Round	Ref. Upper	Ref. Lower
.014"	CE046U	CE046L
.016"	CE047U	CE047L
.018"	CE048U	CE048L

Rectangular	Ref. Upper	Ref. Lower
.014"X .025"	CE051U	CE051L
.016"X .016"	CE049U	CE049L
.016"X .022"	CE050U	CE050L
.017"X .025"	CE052U	CE052L
.019"X .025"	CE053U	CE053L

## Natural Copper Niti Archwire




Super-elastic Niti copper archwire, with thermal properties of wire. Smooth and easy to bend. Once the archwire has been placed, it will return to its original shape and gradually stiffen, aiding your orthodontic treatment as you can alter the force of the archwire in the mouth.


Contain stops

**Pack: 2 pcs**

Round	Natural	
	Without Stops - Ref.	With Stops - Ref.
.012"	CN012	CN012-S
.013"	CN013	CN013-S
.014"	CN014	CN014-S
.016"	CN016	CN016-S
.018"	CN018	CN018-S

Rectangular	Natural	
	Without Stops - Ref.	With Stops - Ref.
.014"X .025"	CN01425	CN01425-S
.016"X .022"	CN01622	CN01622-S
.016"X .025"	CN01625	CN01625-S
.017"X .025"	CN01725	CN01725-S
.018"X .018"	CN01818	CN01818-S
.018"X .025"	CN01825	CN01825-S
.019"X .025"	CN01925	CN01925-S
.020"X .020"	CN02020	CN02020-S
.021"X .025"	CN02125	CN02125-S





## Preformed Archwires

### Natural Dimpled Niti Archwire



Upper  
Lower


- Heat-activation dimpled Niti archwire bent mesially.
- Progressive increase of residual forces with full recovery of shape.
- Specially for self-ligating brackets.

**Pack: 1 pc**

Round	Ref. Upper	Ref. Lower
.012"	DI12U	DI12L
.014"	DI14U	DI14L
.016"	DI16U	DI16L
.018"	DI18U	DI18L

Rectangular	Ref. Upper	Ref. Lower
.016"X .016"	DI1616U	DI1616L
.016"X .022"	DI1622U	DI1622L
.017"X .025"	DI1725U	DI1725L

### TMA Niti Archwire




- TMA - Beta-titanium archwires.
- Superior elasticity compared to steel.
- High surface roughness.
- High friction.
- Ovoid shape.

**Pack: 5 pcs**

	Ovoid	
Round	Ref. Upper	Ref. Lower
.014"	TMA14U	TMA14L
.016"	TMA16U	TMA16L
.018"	TMA18U	TMA16L

	Ovoid	
Rectangular	Ref. Upper	Ref. Lower
.014" x .025"	TMA1425U	TMA1425L
.016"X .022"	TMA1622U	TMA1622L
.017"X .025"	TMA1725U	TMA1725L
.018"X .025"	TMA1825U	TMA1825L
.019"X .025"	TMA1925U	TMA1925L



## Elastic bands



- Latex-free elastics.
- 40 bars of 26 pcs each.

**Pack of 1,040**

Colour	Ref.	Colour	Ref.
Assorted	ME120-C	Orange	ME119
Silver	ME120	Black	ME120-N
Clear	ME113	Green	ME114
Yellow	ME120-A	Red	ME115
Dark blue	ME118	Fuchsia	ME117
Metallic blue	ME122	Light pink	ME116
White	ME120-B	Cyan	ME112
Ivory	ME120-M	Baby Blue	ME112-A
Purple	ME112-M		



- Extra flexible long-lasting elasticity.
- No swelling from moisture in the mouth.
- Normal model 0.24"

**4.57m roll**

	Closed 2.8 mm	Short 3.5 mm	Long 4.0 mm
Ref. Clear	ME130	ME131	ME132
Ref. Grey/silver	ME130G	ME131G	ME132G



Long




Short



Closed







## Elastic bands

### Elastic Separator



- Separators.
- High resistance to breakage.
- Radiopaque.

**Pack of 1,050**


Ref. ME133

### Intraoral elastics




- 100 elastics per bag.
- 50 bags per box.

**Box of 5,000**



Strength	Size				
	1/4" Ref.	1/8" Ref.	3/8" Ref.	3/16" Ref.	5/16" Ref.
Medium 3.5 OZ	EL14	EL18	EL38	EL316	EL516
Medium 4.5 OZ	EL14M	EL18M	EL38M	EL316M	EL516M
High 6.5 OZ	EL14H	EL18H	EL38H	EL316H	EL516H





## Orthodontic Accessories

### Lingual Button



For orthodontic treatments requiring accelerated rotational movements of the tooth or teeth.

Shaped like a button to hook onto the archwire, it is placed in the lingual part of the mouth to perform the rotation with precision and efficiency, providing greater comfort to the patient, ensuring perfect bonding. Secure fixation thanks to its 80-mesh base.

**Pack of 10**

Ref. ME108R Rectangular

Ref. ME108 Round

### Round Aesthetic Lingual Button




Ceramic lingual buttons for clinical cases requiring accelerated rotational movements of the tooth or teeth. Shaped like a button to hook it onto the archwire, it is placed in the lingual part of the mouth to rotate with precision and efficiency, providing greater comfort for the patient, ensuring perfect bonding, as well as great aesthetics.

Rounded base.

**Pack of 10**

Ref. ME108ST





## Crimpable Hook



A very useful accessory in surgical orthodontics to facilitate the handling of the parts during surgery. The hooks are also attached to the bracket to ensure locking of the archwires when placing the elastic bands.

Easy placement in both intraoral and extraoral archwires.

- Secure fixation.
- Maximum arch thickness 0.022"x0.025".

### Pack of 10

	Ref.
Long	ME109
Medium	ME110
Short	ME111

## Extrusion hook



Traction chain for fixed teeth, ideal for inclusion of the affected teeth.

With mesh surface for optimum bonding.

Activation by removal of individual links in the chain.

- Length: 2.5 cm
- Rounded base.

### Pack of 2

Ref. ME208 Gold

Ref. ME207 Silver

## Space Closing Niti Spring



Devices to obtain compressive forces, suitable for closing spaces or to keep them closed temporarily during orthodontic treatment.

Manufactured in Nickel Titanium, they have a great capacity for adaptation, releasing their energy in the form of a light and constant force allowing activation intervals.

- Excellent devices to obtain traction or compressive forces.
- Easy to place and remove.
- Maintain a constant force.

### Pack of 10

	Ref.
0.010 X 6 mm	ME102C
0.010 X 9 mm	ME103C
0.010 X 12 mm	ME104C
0.012 X 6 mm	ME105C
0.010 X 9 mm	ME106C
0.010 X 12 mm	ME107C

## Space Opening Niti Spring




Devices to obtain traction forces, for opening, separating spaces or to keep them open temporarily during orthodontic treatment.

Manufactured in Nickel Titanium, they have a great capacity for adaptation, releasing their energy in the form of a light and constant force allowing activation intervals.

- Excellent devices to obtain traction or compressive forces.
- Easy to place and remove.
- Maintain a constant force.

### Pack of 2

	Ref.
0.010 X 180 mm	ME106
0.012 X 180 mm	ME107



## Orthodontic Accessories

### Archwire Stop



Stop for round and rectangular archwires. Fixation with Weingart pliers.

**Pack of 10**

Ref. ME101

### Ligature Wire Reel



Ligature wires on a reel. For fixing archwires to brackets

**Pack of 4.5m**

	Ref.
0.20 X 30g (0.08mm)	ME127
0.25 X 30g (0.10mm)	ME128
0.30 X 30g (0.12 mm)	ME129

### Class III Adjustable Mask



Ref. MAS3-R (Red)

Ref. MAS3-A (Blue)

- Suitable for children and adults.
- 3-way adjustment.
- Ergonomic.
- Cushioned pads for maximum comfort.
- Comes with replacement pads and Allen key.

**Pack of 1**



## Typodonts

Manufactured in highly durable polymer. Perfect for cementing brackets, for demonstrating to patients. Fixed teeth.

**Pack of 1**

**Typodont without brackets**  
Ref. TI001



**Typodont with metal brackets**  
Ref. TI002



**Typodont with ceramic brackets**  
Ref. TI003



**Typodont with self-ligating brackets**  
Ref. TI004




**Typodont with self-ligating Plus brackets**  
Ref. TI009



**Typodont with sapphire brackets**  
Ref. TI005





## Photography mirrors

Extremely tough scratch-resistant superglass. Made in float glass, specially adapted to take photos, guaranteeing a distortion-free surface. A layer of chrome on both sides delivers maximum-quality reflection. Good surface grip. No duplication of images, 1:1 reproduction.

- Autoclavable up to 220°C, dry and cool air.
- 3 mm thickness.

### Occlusal

**Pack of 1**

	Ref.
Child	ME099
Adult	ME100
Adult XL	ME096



### Lingual

**Pack of 1**

Ref. ME098



### Lateral

**Pack of 1**

Ref. ME097



## Orthodontic Relief Wax



Protective wax to relieve pain and irritation caused by orthodontic appliances (brackets, wires, auxiliary hooks, tubes and ligatures) by sticking to them.

Also recommended for patients who use partial dentures, where wire traction may cause discomfort, sores or abrasions.

Each unit contains 5 orthodontic wax sticks.

**Pack of 50**

Ref. ME214

## Orthodontic Boxes




Orthodontic boxes made of high-quality plastic. For storing orthodontic appliances.

These boxes are very hard-wearing, allowing patients to transport and store their orthodontic appliances safely.

Size: 88 x 82 x 28mm.

**Box of 12 in assorted colours**

Ref. ME215



## Straight Bracket-Removing Pliers



For the removal of all types of bracket. Made in high-quality stainless steel.  
Length 13 cm

Ref. 1537

## TC Distal End Cutter



Especially suitable for orthodontics and prosthetics. Active part in tungsten and small head. For square wire, maximum 0.7mm. For round wire, maximum 0.5mm.  
Length 13cm

Ref. 1538-TC

## Straight TC pliers



Especially suitable for cutting hard wire. Maximum cutting capacity for:

- Square wires 0.41 x 56mm (0.016"x0.022")
- Round wires, maximum 0.7mm (0.020").

Made of German stainless steel with TC inserts.  
Instrument length: 13cm

Ref. 1539-TC

## Straight Cutting Plier 14cm. TC



Specially indicated for cutting fine wire.

Maximum cutting capacity for:

- Round wires 0.012" (0.30mm), 0.014" and 0.016"
- Square wire: 014 x 025" / 016x016" / 016x022".

Made of German stainless steel with TC inserts.

Instrument length: 14cm.

Ref. 1811-TC

## Weingart TC pliers



Specially designed for placing the archwire and for easier bending of the ends. Fine, angled tips to aid access to all areas. Made of stainless steel with TC inserts.

Length: 13cm

Ref. 1540-TC

## TC Micro Ligature Cutters




Ligature cutters. Suitable for orthodontics and prosthetics. Tip with active part made of tungsten, a material that guarantees a long cutting life due to its high resistance to wear. Made in stainless steel. Available with standard or micro tip.

Length 12.5 cm. Length 0.03 mm.

Ref. 1541-TC Standard-tip TC Ligature cutters

Ref. 1400-TC Micro-tip TC Ligature cutters

## Small Circular Punch Aligner Pliers



Piercing pliers to make small circular notches in the exact position. Aligner-cutting pliers to make clean-edged semi-circle cuts in aligners to press the tooth correctly. Also suitable for trimming bothersome areas of gum.

Orthodontic aligner pliers made in high-quality stainless steel.

Length: 13 cm

Ref. 2310

## Circular Punch Aligner Pliers



Perforating pliers to make small crescent-shaped circular notches around buttons.


No need to stick button on the aligner.

These perforating pliers also help eliminate excess plastic so that the soft tissue is not damaged by the splint.

Made in high-quality stainless steel.

Length: 13 cm

Ref. 2311



## Aligner Perforator Pliers



Perforating pliers to exert vertical pressure and achieve a vertical notch in the aligner. Serves as an aid in rotation movements. The vertical pressure on the mesial and distal area of the tooth, both vestibularly and palatally, strengthens the rotation-correcting action of the aligner. Made in high-quality stainless steel.

Length: 13 cm

Ref. 2312

## Keyhole Aligner Pliers



Aligner-perforating Pliers.

These aligner pliers exert pressure and creates a keyhole-shaped notch which is used to anchor extrusion bands, and Class II, Class III, crossed, Z, etc. elastic bands.

Made in high-quality stainless steel.

For optimum maintenance of the pliers, apply a drop of lubricant after sterilisation in the autoclave.

Length: 13 cm

Ref. 2313

## Rectangular Aligner Pliers



Aligner-cutting Pliers

Torque (horizontal) pliers, for creating rectangular notches in the aligner.

These horizontal-pressure pliers aid in root torque movements. Orthodontic aligner pliers made in high-quality stainless steel.

Length: 13 cm

Ref. 2314

## Pliers for bending orthodontic arches 13cm.



Its slim design allows for quick and easy bending of the distal ends inside the patient's mouth.

The angled shape minimizes the friction on the mucosa of the cheeks, which provides relief to the patient and at the same time, allows access to narrower and difficult to reach places, which speeds up and facilitates the work of the professional.

Made of stainless steel with a length of 13cm.

Ref: 2012

## Torque pliers (Male) 13cm.



Specially indicated for torque application, to apply the torque of a rectangular wire in the groove of a bracket and thus alter the inclination of the tooth. To achieve the desired torque, it is necessary to use the two torque pliers (male and female).

For wires up to Ø0.56-0.70mm or 0.022 to 0.028".

Made of stainless steel with a length of 13cm.

Ref: 2013

## Torque pliers (Female) 13cm.



Specially indicated for torque application, to apply the torque of a rectangular wire in the groove of a bracket and thus alter the inclination of the tooth. To achieve the desired torque, it is necessary to use the two torque pliers (male and female).

For wires up to Ø0.56-0.70mm or 0.022 to 0.028".

Made of stainless steel with a length of 13cm.

Ref: 2014

## Band-removing pliers



Orthodontics and prosthetics.


Active parts: one straight part with plastic stop and another, with a short edge that starts from the occlusal pad.

Made in German stainless steel.

Length: 13 cm

Ref. 1535

## SKS Aligner Pliers Kit




Kit contains:

- 1 Circular Punch Aligner Pliers Small
- 1 Circular Punch Aligner Pliers
- 1 Perforating Aligner Pliers
- 1 Keyhole Aligner Pliers
- 1 Rectangular Aligner Pliers


Ref. 23KIT





## Instruments

### Steel Straight Mathieu Needle-holder



14 cm cross-serrated needle holder. Straight Mathieu. Suitable for both surgical and general orthodontics. Length 14 cm, has a spring between the arms and a zip-like fastener so it holds high pressure. Its small grooves on the straight, fine active part also provide improved grip. Made in high-quality stainless steel.

Ref. 465

### Straight TC Mathieu Needle-holder



14 cm cross-serrated needle holder. Straight Mathieu. Suitable for both surgical and general orthodontics. Length 14 cm, has a spring between the arms and a zip-like fastener so it holds high pressure. Its small grooves on the straight, fine active part also provide improved grip. Active part made of tungsten, a material that guarantees a long cutting life due to its high resistance to wear.

Ref. 465TC

### Steel Curved Mathieu Needle-holder



14 cm cross-serrated needle holder. Curved Mathieu. Suitable for both surgical and general orthodontics. Length 14 cm, has a spring between the arms and a zip-like fastener so it holds high pressure. Its gently curved active part has small grooves to improve grip. Made in stainless steel.

Ref. 50

### Examination Kit



Essential kit for a basic dental examination. Stainless steel. Ideal for students.

Kit contains 3 instruments:

- Double probe
- Tweezers
- Cone Socket mirror handle (long thick screw)

Ref. 2000





## Ligature Guide



Guide for ligature placement. Preformed archwire pusher.

An orthodontic instrument specially designed to direct ligatures. Dual active part: one straight and one curved, for maximum access.

Made in high-quality German stainless steel.

Ref. 2110

## Curved tweezers for Tubes



Tweezers for placing buccal tubes in Orthodontics. Crossed, curved active part ending in a very fine tip. Stainless steel.

Ref. 2001

## Tweezers for Brackets



Tweezers for placing brackets.

Short and flat active part for placing brackets. Stainless steel.

Ref. 401

## Star-shaped Bracket Placer



Universal bracket placer.

To measure the height to place the slot. 4 heights (3.5, 4, 4.5 and 5mm). Stainless steel.

Ref. 397

## Bracket Position Measurer/Placer



Measurer for bracket positioning.


To measure the correct distance between the bracket and the incisal edge.

Dual active part ending in swivel head.

For heights of: 4.0 - 4.5 and 3.5 - 5.0 mm. Stainless steel.

Ref. 1424





## Ligatures and Wires

### Kobayashi Long Metal Ligatures for Orthodontics



- Kobayashi ligatures in high-quality stainless steel.
- For fast placement with traction hook in hookless brackets.

**Pack of 100**

	Ref.
0,012	ME206
0,014	ME209

### Kobayashi Short Metal Ligatures for Orthodontics



- Kobayashi ligatures in high-quality stainless steel.
- With blind hook point and screw tip.

**Pack of 100**

	Ref.
0,012	ME204
0,014	ME205

### Preformed Long Metal Ligatures for Orthodontics



- Preformed ligatures in high-quality stainless steel.

**Pack of 100**

	Ref.
0,010	ME210
0,012	ME211

## Preformed Short Metal Ligature for Orthodontics



- Preformed ligatures in high-quality stainless steel.
- With long braid for use with ligature instruments.

**Pack of 100**

Ref. ME213

## Retention Wire .028" Long



Long retention wire.

- Length: 0.014"

**Pack of 20**

Ref. ME212

## 6 Thread Retention Wire



Passive 6-thread retention wire 8 pcs. 6-thread braided wire rod with for lingual retention. For manufacturing permanent fixed lingual retainers.

**Pack of 10**

Ref. ME216





## 2.5L Ultrasonic Bath



Recommended for cleaning small items in dental laboratories and clinics. Excellent cleaning of removable prostheses.

Simple and user-friendly. Heating system with electronic timer.

The removable tank makes cleaning easier and improves hygiene.

Timer setting.

WARRANTY: 1 year

Ref. AS004



## 2.5L Ultrasonic Bath



Medium ultrasonic bath.

Convenient and practical for cleaning small items in dental laboratories and clinics.

Its 2.5-litre tank has the right dimensions to hold the usual instruments used in the dental clinic.

With heating system with electronic timer.

Straightforward and easy to use, achieves excellent cleaning of removable dentures using ultrasonic-frequency waves.

The removable tank makes cleaning easier and improves hygiene. Supplied with a protective plastic tray for the tank.

With timer setting.

WARRANTY: 1 year

Ref. AS004-N





## Distillers / Gutta Cutter

### Water Distiller Jug



#### DIMENSIONS:

- Length: 23.5cm
- Width: 22cm
- Height: 38cm

Capacity to produce 4 litres of water in approximately 4 hours (depending on the room temperature), which means a daily average of 20 litres of water. Ideal for domestic and corporate use.

WARRANTY: 1 year

Ref. AS010

#### TECHNICAL SPECIFICATIONS:

- Stainless steel water tank.
- Comes with reusable spare filter.
- Nominal voltage: 220-240V - 50-60Hz
- Nominal power: 750 W
- Capacity: 4L
- Distillation volume: 1.5L / hour
- Protection from overheating: shuts down the appliance

### C-Blade Gutta Cutter



Facilitates the cutting of gutta-percha tips at the entry to the endodontic canal. The end of the Gutta Cut is inserted around the excess gutta percha and by pressing the button, the tip heats up, producing a clean, fast, precise and safe cut.

WARRANTY: 1 year

Ref. C-BLADE


#### FEATURES:

- Ergonomic design of the cutting end with four interchangeable tips of various shapes, for wide or narrow canals.
- Autoclavable tips.
- Do not immerse in disinfectant liquids.
- Lightweight and wireless.
- 80 cuts per charge.

#### CONTENTS:

- One Gutta Cut gutta percha cutter
- 1 Charger
- 1 rechargeable battery
- 4 tips





## Electric Heat-Sealing Machine



Easy to use, compact and ergonomic. It has an LCD screen to display procedures being performed at all times.

WARRANTY: 2 years.

Ref. AS011E

### TECHNICAL SPECIFICATIONS:

- Standard Dimensions: 539 x 517 x 205mm
- Dimensions with double roll holder: 539 x 517 x 205mm
- Supply voltage: 200-240 V; 50Hz.
- Sealing width: greater than/equal to 12mm
- Maximum sealing length: 300mm
- Weight: 8kg.

LCD screen  
to view procedures



[www.sksdental.es](http://www.sksdental.es)

*Our Values  
serving you online*



**MAXIMUM  
SAVINGS**



**DELIVERY  
WITHIN 24-48 HOURS**



**UNBEATABLE  
DEALS**



**FREE SHIPPING**  
for purchases over €100  
(excl. VAT)




**SAFE SHOPPING**



**15 DAYS FOR  
REFUNDS**  
(except for equipment)







## Contra-angle Handpieces and Turbines

\* Rotary instruments from AS Technology offer a 1-year warranty on the part and 6 months on the rotor.

### Turbine with light + Multiflex Attachment\*



Speed of up to 310,000 rpm. Despite its power, it is an auxiliary turbine and not suitable for intensive use. Its LED lighting enhances visibility of the area to be treated. Large head and 360° rotation making it more comfortable to use by significantly reducing fatigue during treatment. With triple spray cooling. Push button system for quick and effortless changing of burs. Autoclavable at 135°C.

Ref. AS012: Turbine sheath with light + multiflex attachment

Ref. AS008: As Technology FK600 Turbine with light

Ref. AS005: As Technology FK500 Multiflex attachment

#### FEATURES:

- Auxiliary turbine with light.
- Not suitable for intensive use.
- MULTIFLEX.
- LED lighting
- 360° rotation.
- Speed of up to 310,000 rpm.
- Triple spray.
- Multiflex Light connection.
- Bur-changer button.
- Autoclavable at 135°C.

### 1:5 Multiplying Contra-angle with light\*



1:5 multiplying CA rotary instrument with internal spray. Suitable for implantology. Change burs by pushing button. Optimum speed of up to 200,000 rpm. Convenient and practical during treatment, its ergonomic design makes it efficient. Easy to clean. Autoclavable at 135°C. With light.

Ref. AS007

#### FEATURES

- Multiplier contra-angle.
- 1:5 Multiplier CA.
- For implantology.
- Internal spray.
- Maximum speed up to: 200,000 rpm.
- Convenient, practical and efficient.
- Bur-changing button.
- Autoclavable.
- With light.

### 1:1 Contra-angle with light\*



High-quality contra-angle handpiece for micromotor with steel body and handle and interchangeable head, incorporating a bur-changer button. Autoclavable at 135°C.


1:1 with internal cooling. Made in stainless steel. Speed of up to 25,000 rpm.

Ref. AS001

#### FEATURES

- Max. speed 25,000 rpm.
- Internal spray.
- Sterilisable and thermo-disinfectable.
- 1:1 transmission.
- With bur-changer button.
- Ergonomic and hygienic design.





## 2g syringe refills.

Colour	Ref.
A1	004-101A1
A2	004-101A2
A3	004-101A3
A3.5	004-101A3.5

## Flow Composite Kit 4 x 2g syringes and refills

Low-viscosity light-curing fluid resin, highly aesthetic, used for all types of cavity preparations, and with excellent polishability.

Easy to polish, high-gloss, high resistance to compression, wear- and stain-resistant. Radiopaque to facilitate identification on x-rays.

The composite has a particle size of approx. 0.7 microns. An ideal option for use in the following cases:

- Small anterior and posterior restorations.
- As a pit and fissure sealant.
- Small defects and filling of cervical gaps caused by erosion from brushing or gum recession.
- Eliminate retentions before taking inlay impressions.

### ADVANTAGES AND FEATURES:

- High-strength micro-hybrid formula with excellent polishability.
- Fast, light-curable in just 30 seconds for light shades; 40 seconds for dark shades.
- Radiopaque
- Fluoride leaching.
- Low contraction.
- High resistance to compression.
- 2g syringes.

### Kit contents:

- 4 x 2g syringes. A2, A3, A3.5, B2
- 20 dispensing tips

Ref. 004-010

## Micro-hybrid Kit of 4 x 4.5g syringes and refills



## 4.5 g syringe refills.

Colour	Ref.
A1	001-401A1
A2	001-401A2
A3	001-401A3
A3.5	001-401A3.5

Ideal composite for anterior and posterior restorations.

Inorganic filler.

0.7 micron average particle size


### ADVANTAGES AND FEATURES:

- Easy to manipulate during placement.
- Does not stick to instruments.
- Excellent polishing properties.
- High aesthetic result and high gloss. Ideal composite for anterior and posterior restorations.

### Kit contents:

- 4 x 4.5g syringe. A2, A3, A3.5, B2
- 7ml bond.
- 7ml surface sealant
- 4g acid etch.
- Accessories.

Ref. 001-414



## Acid etch

### Prime-Dent Etching Gel/Fluid

37% phosphoric acid, for etching during bonding procedures.

Provides even etching, creating the ideal surface for perfect bonding of resins.

Prime-Dent Etching Gel is thixotropic, stays in place and does not run.

### Etching gel kit 4 x 1.2ml syringes.

**Set of 4 syringes.**

Ref. 008-010SF



### 50 ml Etching Gel Syringe

**Pack contains:**

- 50 ml syringe
- 50 dispensing tips
- 3 empty syringes with attachments

Ref. 008-060SF




### 14g Etching Gel Syringe

**Pack contains:**

- 14g syringe.
- 3 dispensing tips

Ref. 008-030SF





## Dual Cure Dentine/Enamel Bonding Resin 7ml.



Fifth-generation multi-purpose, hydrophilic light-curing bond, suitable for use on moist dentine (wet-bonding technique).

The adhesive incorporates bonding elements all in one tub, designed to bond composites and compomers to dentine and enamel, as well as treated amalgam, metal and ceramics.

After etching, penetration and saturation of dentine and enamel is carried out in 3 layers.

**Pack containing 7ml.**

Ref. 006-030

## Light Cure Self Etch Prime & Bond 7ml.



Light Cure Self is a state-of-the-art adhesive.

No mixing needed as it incorporates etching and bonding elements in the same container.

Quick and time-saving. High quality, reduces sensitivity and increases bond strength.

**Pack containing 7ml.**

Ref. 006-033

## Dual Cure Luting Cement 10g Syringe



### Kit contains:

- 1 x 10g syringe.
- 5 mixer tips
- 5 intraoral tips


Ref. 100-101

Dual cement for use in final cementing of implant-supported prostheses, crowns, bridges, inlays, onlays and veneers. Cementing of posts, screws and periodontal splinting.

Glass ionomer composition in a Bis-GMA matrix base of dental resins, activators, catalysts and pigments.

### Advantages and Features:

- Rapid cementing.
- Auto-Mix syringe.
- Practical, economical and especially useful for multiple surgical procedures.
- Free from air bubbles, no wasted material due to surplus mix.
- Predictable, accurate, consistent and homogeneous mixing at all times.
- Minimal heat generated.
- Protects the pulp.
- Radiopaque.
- Easy to identify on x-rays.
- Low thickness film.
- Provides easy crown and bridge fitting
- Fluoride leaching.
- Reduces secondary caries.



### 10g syringe Ionomer band cement

Dual resin-and glass-ionomer-based cement.

Ideal for cementing metal bands as it chemically bonds the tooth structure to the metal band creating a strong anchorage.

Excellent consistency.

Dual setting with fluoride leaching.

Quick cementing, easy to clean and ideal for bite raising.

Polymerisable to speed up the process.

No need for light-curing.

Full chemical process in five minutes.

#### Kit contains:

- 1 x 10g syringe.
- 5 mixer tips
- 5 intraoral tips
- Instructions.

Colour	Ref.
Blue	015-061
Natural	015-028



### Light Cure Kit Orthodontic Bond 1x5g.

Light-curing, single-component resin system for bonding all types of supports to enamel.

Especially suitable for bonding metal, resin or ceramic brackets.

Precise viscosity to prevent bracket slippage.


#### Kit contains:

- 1 x 5g syringe of paste.
- 1 x 3ml bond.
- 1 x 3ml etch gel.

Ref. 012-022







## Non-sterile Cotton Rolls



Cylinders for dental use, soft and absorbent

Composition: 100% cotton

Store in a dry and ventilated place.

**Packaging with 300 g.**

Size	Ref.
N°1 - 8x38 mm.	ME1301
N°2 - 10x38 mm.	ME1302
N°3 - 12x38 mm.	ME1303

## Disposable Swabs

Disposable swabs to apply the adhesive or other materials such as sealants, flowable composites, varnishes on the dental surface.

**The pack contains 100 swabs.**





Thickness	Color	Ref.
Fine	<span style="color: yellow;">■</span> Amarillo	SKS002
Regular	<span style="color: blue;">■</span> Azul	SKS003




kike  
TOYS  
by SKS®

When going to the dentist  
becomes child's play



The world  
of Kike



# Fantomas



Fun, cute cuddly toys, specially designed for little ones to help explain how to brush and clean their teeth properly.

- ◆ Can be used with dental floss.
- ◆ Come with large toothbrush.
- ◆ Ideal for cementing brackets onto the toy.
- ◆ Suitable for use with 3-in-1 syringe and water in their mouth.



Pink Panther

REF: KT\_039



Hippo

REF: KT\_047



## Dragon

REF: KT\_040



◆ Glove puppet.



## Lion


REF: KT\_041





## Monkey

REF: KT\_043



## Dinosaur

REF: KT\_042



◆ Glove puppet.



## Leopard

REF: KT\_044



## Molar Cushion

Cuddly cushion in the shape of a molar.

Ideal for creating the right atmosphere for your clinic and entertaining children in your waiting room.

Cuddly toy toothbrush

REF: KT\_046

Cuddly toothpaste cushion

REF: KT\_045





# Backpacks




## Butterfly

Neoprene fabric

Dimensions: 33 x 26 x 12cm



Yellow  
REF: KT\_002





### Bee

Neoprene fabric

Dimensions: 21.5 x 21cm


REF: KT\_007

### Dinosaur

Neoprene fabric

Dimensions: 30 x 23.5 x 14cm

REF: KT\_004





### Shark

Dimensions:  
28.5 x 24.5 x 9 cm

REF: KT\_005




### Frog

Neoprene fabric

Dimensions:  
21.5 x 21 cm

REF: KT\_006





# Toothbrushes

## Children's toothbrushes

- ◆ Extra-soft bristles to be gentle on children's gums.
- ◆ Ergonomic handle
- ◆ Recommended for 6 years and up



Box: 12 pcs, assorted  
REF: KT\_033



# Off to school

## Kike Toys Pencils

Pencils with fun designs of all our friends from Kike Toys. With eraser tip.

### Bag of 48

- Kike Toys Blue Little Friends Pencils  
REF: KT\_017
- Kike Toys Pink Little Friends Pencils  
REF: KT\_018



## Super Patient Pencils

Kike Toys Draki design. With eraser tip.

### Bag of 48

REF: KT\_019





## School pencil cases

Pencils with fun designs of all our little friends from Kike Toys.  
With transparent pocket to put clinic card in and personalise it.

Sold individually or with contents.

### Pencil cases with contents.

#### Contains:

- 2 molar-shaped erasers
  - 1 Kike Toys design pencil
  - 1 Children's toothbrush
  - 1 Children's hourglass
- Dragon design pencil case with contents  
REF: KT\_023
  - Giraffe design pencil case with contents  
REF: KT\_024



### Pencil cases without contents

#### Pack of 25 without contents:

- Dragon design pencil case  
REF: KT\_025



- Giraffe design pencil case  
REF: KT\_032





## Molar erasers

Fun molar and smiley face shaped erasers from Kike Toys.

- Bag of 100 white molars  
REF: KT\_020
- Bag of 200 molars, assorted colours  
REF: KT\_021

## Smiley erasers

Fun molar and smiley face shaped erasers from Kike Toys.  
Pack of 50

- Blue smiley eraser  
REF: KT\_050
- Green smiley eraser  
REF: KT\_051
- Yellow smiley eraser  
REF: KT\_052
- Pink smiley eraser  
REF: KT\_053



## Assorted Best Patient Certificates

Best patient certificates designed with all of our little friends from Kike Toys.

A practical gift for our bravest and keenest little ones.

Contains 20 pcs, assorted  
(10 of each model).

REF: KT\_038





# Playtime

## Assorted stickers

Assortment of stickers of our little friends from Kike Toys.

Contains 100 pcs, assorted

REF: KT\_030




## Smiley Face Ink Stamps

Ink stamps with various designs of funny faces.

Ideal distraction for children in the surgery.

Pack of 24, assorted

REF: KT\_057



# Playtime

## Tooth-holder Chests

In the shape of a treasure chest, for storing milk teeth.

Assorted Colours

Made from plastic.

Dimensions of each chest: 4 x 2 x 2cm

Bag of 50

REF: KT\_055





## Tooth-holder Mice

Fun little mice for patients to keep their milk teeth in.  
Black flexible rubber ears and tail.

Bag of 50 assorted

REF: KT\_013



## Paste and Stress-relief molars Kit

Fun, relaxing molar- and toothpaste-shaped polyurethane sponge toys.

Help to reduce stress and anxiety levels by squeezing and releasing the toy in the hand.

Ideal as a gift for children to entertain paediatric patients during dental appointments.

Stress-relief molars with different screen-printed faces.

Dimensions:

Toothpaste:

15cm (high) x 4cm (wide)

Molars: 1 cm (high) x 4cm (wide)

Contains 4 pcs:

- 1 blue toothpaste
- 1 pink toothpaste
- 1 blue molar
- 1 pink molar

REF: KT\_058





# Playtime

## Tooth Wristbands

Cord wristbands with smiley teeth in assorted colours.

Bag of 48 assorted

REF: KT\_011



## Hourglasses

Smiley hourglasses in fun colours.

Bag of 10

- ◆ Yellow REF: KT\_016
- ◆ Red REF: KT\_022
- ◆ Purple REF: KT\_028
- ◆ Blue REF: KT\_012
- ◆ Green REF: KT\_031
- ◆ Orange REF: KT\_049



## Rainbow Slinky Springs

Mini rainbow-coloured spiral springs.

Closed spring dimensions:

Diameter: 3 cm

Height: 3.3cm

Packaging dimensions:

19cm wide x 14 cm deep x 3cm high

Not recommended for children under 3 years.

■ **Box of 24**

REF: KT\_010

## Funny Face Bouncy Balls

Funny face bouncy balls

Assorted colours and shapes, shiny and bright.

Durable plastic. Ø 3cm.

Not recommended for children under 3 years.

■ **Bag of 48**

REF. KT\_027



# Mirror, mirror...

## Orthodontic Dental Hand Mirror

Kike Toys Molar-shaped mirror in the shape of a tooth, made from methacrylate, with a fun design for orthodontics.

Good fit to the hand and broad viewing surface.

Dimensions: 14cm wide x 24.5cm high

■ Pink REF: KT\_015

■ Blue REF: KT\_014



Be bracket  
my friend




## Dental Hand Mirror

Kike Toys Molar-shaped mirror made from methacrylate.

Good fit to the hand and broad viewing surface.

Dimensions: 14cm wide x 24.5cm high

- ◆ Red REF: KT\_034
- ◆ Blue REF: KT\_036
- ◆ Green REF: KT\_037
- ◆ Yellow REF: KT\_035
- ◆ Pink REF: KT\_029
- ◆ Purple REF: KT\_054




# SKS Xpress Light Wireless Lamp

**SKS Xpress Wireless Light.**  
**Instant light-curing Lamp**

Ref. 200-009 ORT (Black)

Ref. 200-010 ORT (White)



New **3''**   
seconds  
**curing**  
**Orthodontic Mode**



962 425 366



[www.sksdental.es](http://www.sksdental.es)



[info@sksdental.es](mailto:info@sksdental.es)



# Light-curing lamp

## SKS LIGHT - LED COLOUR CURING LIGHT

### SKS LIGHT - LED COLOR CURING LIGHT Mini LED light-curing lamp.

Contents: 1 Charging adapter with micro USB 5V DC / 1.2 A output, 1 Wireless probe, 1 Battery module, Light guide, 1 Protective filter, 1 Protective cover, 1 Base, Instruction manual.

Ref. Reference 200-003PINK (Pink)

Ref. Reference 200-003BLUE (Blue)

Light wavelength: 410–490 nm

Weight 100 g. Lamp  
276 g. Metal base

Maximum power: 1100 mW/cm<sup>2</sup>,  
Timing: 10 sec, 20 sec, 30 sec.

360° rotation and automatic  
compensation circuit for  
constant power output

Wireless probe dimensions:  
(18 mm minimum diameter  
and 25 mm maximum,  
length without light guide  
148 mm)



Warning!!! ⚠

Products marketed by SARKIMPISA IMPORT S.L.U comply with legal healthcare regulations. SKS DENTAL shall not be liable for any damages resulting from the use of these products, which are intended strictly for professional dental use.

Unsuitable toys and toys that may be dangerous for children under 3 years of age must be used under adult supervision. This catalogue contains items that may include parts that are dangerous for minors.

SARKIMPISA IMPORT S.L.U shall not be liable for any errors in publication, and reserves the right to correct them as appropriate.

Your purchase data, which is necessary to create your order, is collected on a voluntary basis for strictly commercial purposes, in a file managed and received by SARKIMPISA IMPORT S.L.U in accordance with the provisions of Spanish Organic Law 3/2018 of 5 December, on the Protection of Personal Data and Guarantee of Digital Rights, you may exercise your rights of access, rectification, erasure, restriction of processing, portability and objection in accordance with the provisions of Articles 12 to 18 of the aforementioned law, by contacting our customer service department or by sending an email to [info@sksdental.es](mailto:info@sksdental.es).

