



# CONVERSIONOR DE BANDAS



962 425 366



[www.sksdental.es](http://www.sksdental.es)



[info@sksdental.es](mailto:info@sksdental.es)

## PRIMER MOLAR SUPERIOR (6)

SKS DENTAL	UNITEK 3M	GAC	LEONE	ORMCO	LANCER	ORTHO ORGAN	DENTAURUM
29	29						1,2
29+	29+	1		1,2	1,2	1	3,4
30	30	2		3	3,4,5	2	5
30+	30+	3		4,5,6	3,4,5	3,4	6,7
31	31	4		4,5,6	6	3,4	8
31+	31+	5		7	7	3,4	8
32	32	6	8	8,9,10	8	5,6	9,10
32+	32+	7	8	8,9,10	9	5,6	9,10
33	33	8	9	11	10	5,6	11
33+	33+	8		11	11	7	12
34	34	9	10	12,13	12	8, 9	13,14
34+	34+	10	11	14,15	13	10	15
35	35	11,12	12	14,15	14	11,12	16
35+	35+	13	13	17	15,16	11,12	17,18
36	36	14,15, 16	14	17	17	17,18	19
36+	36+	14,15,16	15,16	18,19	18,19	17,18	19, 20
37	37	17	17	20	20	15	21
37+	37+	18	18,19	21	21	16,17	22
38	38	19,20,21	20,21	22,23	22,23	18,19,20	23,24
38+	38+	22	22,23	24	24,25,26	21, 22	25,26
39	39	22	24	24,25	24,25,26	21, 22,23	25,26
39+	39+	23	25	27	27	24	27
40	40	23	26	28	28	25	28,29,30
40+	40+	24	27	29,3	29	27,28,29,30	28,29,30
41	41	25,26	28 Y 29	31	31	31	31
41+	41+	27,28,29,30	30	32	32	32	32
42	42	0	31	32	30		
42+	42+	0	0	0	31		
43	43	0	0	0	32		



## PRIMER MOLAR INFERIOR (6)

SKS DENTAL	UNITEK 3M	GAC	LEONE	ORMCO	LANCER	ORTHO ORGAN	DENTAURUM
29	29		1				
29+	29+	3,4,5	2,3	1			1
30	30	6,7	4	2	1,2	1	2
30+	30+	8	5	3,4	3,4	2	5
31	31	9	6	5,6,7	5	3,4	6
31+	31+	10	7	5,6,7	5	3,4	
32	32	11	8	8	7,8	5,6	9
32+	32+	11	8	9	9	5,6	10
33	33	12	9	10,11	10	7	12
33+	33+	13	9	13	11	8	13
34	34	14	10	14	12	9	14
34+	34+	15	11	14	13	9	15
35	35	16,17	12	16,17	14	10	16
35+	35+	18	13	16,17	15,16	11	17,18
36	36	19,20	14	19	17	12	19
36+	36+	21,22	15,16	20	18,19	13,14	20,21
37	37	23	17	21	20	15	22
37+	37+	24	18,19	22,23	21	16	23
38	38	25	20,21	24	22,23	17,18	24
38+	38+	25	22,23	25	24,25,26	17,18	25
39	39	26	24	26	24,25,26	19	26,27
39+	39+	27	25	27	27	20	28
40	40	28,29	26	28,29	28	21,22,23	29
40+	40+	30,31	27	28,29	29	24,25,26	30
41	41	32	28 Y 29	30	31	27,28,29	31
41+	41+	32	30	30	31		32
42	42		31				
42+	42+						
43	43						